



Philips PageWriterTouch



Speed and **accuracy**
set this cardiograph apart

PHILIPS

See and touch the quality

See the Philips PageWriterTouch and get a clear sense of the future of cardiography — where precision performance and sophisticated ergonomics converge.

Prepare to be wowed at your first encounter and for a similar sensation every time you use it. And PageWriterTouch is designed for endless use — in large hospitals where ECG volume is high and speed and quality are paramount.

Simply touch the large (15-in.), high-resolution color display and you're in business — recording, selecting and printing the highest quality ECGs faster than you can imagine. And PageWriterTouch stores ECGs in XML format, so you can transmit to any system that understands this industry standard.

From accurate and complete data entry to comprehensive data fields that prompt users to collect required information, every detail of PageWriterTouch is designed to streamline ECG testing.

And when you're responsible for hundreds of ECGs across miles of hospital terrain every day, attention to detail makes all the difference.



FAST START: The always-ready patient module neatly cradles lead wires, so there's no tangled mess to slow you down.

With PageWriterTouch, expect your every move to be faster and more precise based on its:

- High-resolution, 15-in. color touch display
- Choice of data entry methods
- Comprehensive analysis
- Standards-based operation

And PageWriterTouch is also designed to make every step easier to take:

- **Simple to connect.** Patient module ensures lead wires are tangle-free and ready to use.
- **Simple to operate.** Ergonomic design means any user can comfortably touch the display, type on the keyboard or push the cart.
- **Simple to transmit.** Move data from the cardiograph to your ECG Management System quickly and easily.

Sharpen your senses

Through an ingenious union of form and function, the Philips PageWriterTouch delivers a whole new look, feel and standard for ECG recording.

The only thing not sharp about it are its corners.

PRINT the complete 12-lead ECG, RHYTHM and FULL-DISCLOSURE.

MANEUVER smoothly and comfortably with glide-ready wheels and a handle positioned at just the right height.

SEE AND TOUCH what you need quickly and clearly: lead quality, ECG snapshots, the complete, print-ready ECG and interpretation.

CONNECT via LAN, wireless LAN or modem to your ECG Management System.

ENTER DATA by typing on the real keyboard, or choose the optional magnetic card reader or bar code scanner.

LEAD WIRES are always organized and ready to apply when you lift the patient module.

ORGANIZE and safeguard hospital assets and personal items with one adjustable shelf or more (optional) and the optional locking drawer.



GET TO IT touching the you need to

Raise your expectations



The depth is optimal for the display and one touch is all record.



JUST RIGHT: The height is optimal for pushing and keyboard typing.

This fast and accurate cardiograph, with its advanced interpretive algorithm, delivers for technicians, physicians, patients — and the bottom line:

- **More tests in less time with less effort**
- **Improved documentation**

From the moment you turn it on, the Philips PageWriterTouch contributes to your effectiveness.

Everything inside PageWriterTouch is designed to make ECG acquisition, analysis, printing, storage and transmission a faster and more accurate process.

At its core is the Philips 12-lead Algorithm, a sophisticated program that ensures highly sensitive and specific interpretation.

The culmination of extensive research and development, this is the same algorithm found in the Philips IntelliVue Patient Monitoring System. So far it has passed some of the toughest standards for accuracy and reliability — ours.



The Philips 12-lead Algorithm ensures highly sensitive and specific interpretation.



From full-disclosure review, select any 10 seconds to get a complete, interpreted 12-lead ECG.

PageWriter **touch**


Experience the speed



FAST CONNECT: Quickly and easily attach lead wires from your patient's left or right side. A full-color lead map pinpoints problem electrodes on the torso.



FAST ENTRY: Enter patient data quickly and accurately as you choose: type on the real keyboard, scan the bar code or use the magnetic card reader. You can also download patient data and orders from a floppy disk.



FAST TOUCH: Create a portfolio of ECG snapshots by tapping the screen. Select and touch a thumbnail to preview in full size. Choose the best one to print. Do it all in seconds.

For the most widely used diagnostic test in the world, the standout choice is the Philips PageWriterTouch. We say that with confidence because we developed it based on extensive market input.


And then we delivered more.

From the instant-on feature to the quick connection to your ECG Management System, PageWriterTouch accelerates ECG testing from start to finish.


Speed and accuracy set this cardiograph apart. See for yourself.

Training at your fingertips

To get your entire department up and running on PageWriterTouch, Philips puts training at your fingertips.



A comprehensive computer-based training program ships with every PageWriterTouch. This self-paced program with a pre- and post-test qualifies users for CEUs.



No need to post notes on how to operate the cardiograph. This simple-to-understand screen gets you started in a flash. If you choose, you can quickly and easily replace this start-up screen with one of your own.

Philips Medical Products deserve Philips Supplies

To maximize the performance of your Philips medical equipment and protect your investment, don't settle for anything less than Philips Medical Supplies. Quality and reliability are 100 percent guaranteed. All Philips medical supplies and accessories have passed Philips rigorous safety and performance tests. In addition, we offer full product and technical support. Check out our complete interactive supplies catalog online at <http://shop.medical.philips.com>.

Philips Medical Systems is part of Royal Philips Electronics

INTERESTED?

Would you like to know more about our imaginative products? Please do not hesitate to contact us. We would be happy to provide specific information about our products and services, or put you on our mailing list for news about new product developments, upcoming events or for our clinical journal, MedicaMundi. We would be glad to hear from you.

On the web

www.medical.philips.com

Via e-mail

medical@philips.com

By fax

+31 40 27 64 887

By postal service

Philips Medical Systems

Global Information Center

I.B.R.S. / C.C.R.I. Numéro 11088

5600 VC Eindhoven

Pays-Bas / The Netherlands

(no stamp required)



**Philips Medical
Systems is part
of Royal Philips
Electronics**

www.medical.philips.com
medical@philips.com
fax: +31 40 27 64 887

United States
(800) 934-7372

Canada
(800) 291-6743

Europe, Middle East and
Africa
(+49) 7031 463 2121

Latin America
(954) 835-2600

Asia Pacific
(852) 2821-5888

© Koninklijke Philips
Electronics N.V. 2003
All rights are reserved.

Reproduction in whole or in
part is prohibited without the
prior written consent of the
copyright holder.

Philips Medical Systems
Nederland B.V. reserves the right
to make changes in specifications
and/or discontinue any product at
any time without notice or
obligation and will not be liable
for any consequences resulting
from the use of this publication.

Printed in the United States
July 2003
4522 982 91111