

Holter Compatible Operating Systems

Operating System	Holter 2.9.4	Holter 3.0
Windows XP Professional, SP2	✓	✓
Windows XP Professional, SP3	✓	✓
Windows Vista Business or Ultimate, SP1	✓	✗
Windows 7 Professional (32 Bit), SP1	✓	✓
Windows 7 Professional (64 Bit), SP1	✗	✓
Windows 7 Enterprise (32 Bit), SP1	✗	✓
Windows 7 Enterprise (64 Bit), SP1	✗	✓
Windows Server 2003 Standard Edition, SP2 (32 Bit)	✓	✓
Windows Server 2003 Standard Edition, SP2 (64 Bit)	✗	✓
Windows Server 2008 Standard Edition, SP2 (32 Bit)	✗	✓
Windows Server 2008 Standard Edition, SP2 (64 Bit)	✗	✓

✗ Not supported

Recorder Compatible Operating Systems

Operating System	DT+		DTXT	
	Holter 2.9.4.x	Holter 3.0	Holter 2.9.4.x	Holter 3.0
Windows XP Professional, SP2	✓	✗	✓	✗
Windows XP Professional, SP3	✓	✓*	✓	✓
Windows Vista Business or Ultimate, SP1	✓	✗	✓	✗
Windows 7 Professional (32 Bit), SP1	✓*	✓*	✓	✓
Windows 7 Professional (64 Bit), SP1	✗	✗	✗	✓
Windows 7 Enterprise (32 Bit), SP1	✗	✓*	✗	✓
Windows 7 Enterprise (64 Bit), SP1	✗	✗	✗	✓

✓* = Requires DigiTrak Plus (DT+) Firmware 4.1.11

- DigiTrak Plus (DT+) not supported on 64bit Windows
- DigiTrak Plus (DT+) EOL: 31 Dec 2014
- DigiTrak XT (DTXT) Firmware: C.2